

# CEST MRI Using Golden-Angle Cartesian Acquisition with Compressed Sensing Reconstruction

Ding Xia<sup>1</sup>, Rodolphe Leforestier<sup>1</sup>, Li Feng<sup>1,2</sup>, Xiang Xu<sup>1</sup>

1. BioMedical Engineering and Imaging Institute, Icahn School of Medicine at Mount Sinai, New York, United States
2. Center for Advanced Imaging Innovation and Research (CAI2R), New York University Grossman School of Medicine, New York, United States

## INTRODUCTION:

We propose a new 3D dynamic CEST MRI technique which involves saturation transfer preparation, a gold-angle rotated, variable-density signal acquisition and compressed sensing (CS) reconstruction. The approach offers improved signal detection sensitivity and motion artifacts reduction, enables a continuous signal acquisition and retrospective reconstruction with flexible temporal resolution.

## METHODS:

The sequence framework combines golden-angle rotated, variable-density spiral-centric Cartesian acquisition<sup>1, 2</sup> and low-rank based multi-coil CS reconstruction.<sup>3</sup> Following the preparation pulses, one shot, which includes a number of phase-encoding steps predefined in the ky-kz Cartesian grid, is acquired. Following shots are rotated by a golden angle to allow for a uniform and continuous coverage of k-space. 400 k-space lines were acquired per shot after the CEST preparation within an imaging matrix of 96x96x20. For each frequency offsets, 5 shots were acquired corresponding to the time required for a fully-sampled k-space. The CEST preparation uses a train of 20 Gaussian pulses, 50 ms each, with inter-pulse delay of 0.5 ms, and B1 power was 1.6  $\mu$ T. The sampling trajectory was rotated by 137.5° between different shots. For image reconstruction, first, images were reconstructed by CS reconstruction with spatiotemporal constraints using all 5 shots (Fig. 1a). Then “shots” were retrospectively removed from the k-space and CS reconstruction was applied to the first 2 shots (Fig. 1b) and the first shot alone (Fig. 1c). The PLOF<sup>4</sup> method was used for extracting the NOE contrast at -3.5 ppm from a posterior white matter ROI and generating the NOE maps.

## RESULTS:

The CS reconstruction was able to generate artifact-free images with just 2 shots or a single shot, which correspond to acceleration factors of 4.8 or 9.6 and a temporal resolution of 10s or 5s per 3D volume per frequency offset. The Z-spectral features from multiple exchange pathways were well preserved (Fig. 2). Fitting results showed almost identical NOE contrast from images reconstructed from 5, 2 and even a single shot (Fig.3). The NOE maps showed slight blurring at the tissue boundaries when using a single shot for reconstruction (Fig.3).

## DISCUSSION:

Combining the advantages of our variable density spiral-centric Cartesian acquisition and CS reconstruction, we demonstrated that the CEST contrast can be well preserved with 4.8-fold acceleration rate. The choice of further acceleration depends on the application. For dynamic imaging, e.g., repeated observation at a certain frequency offset, the current sequence framework offers flexibility of retrospectively combining different number of shots to achieve arbitrary temporal resolution.

## CONCLUSION:

We proposed a 3D CEST imaging method using a golden angle rotated and variable-density Cartesian acquisition in combination with multi-coil compressed sensing reconstruction. It is possible to obtain whole brain CEST maps well under 5 min with this method, making it applicable to clinical applications.

## ACKNOWLEDGMENTS:

NIH funding support, R00EB026312, R21AG083685

## REFERENCES:

- [1] Cheng JY, et al. *J Magn Reson Imaging*. **2015**, 42 (2), 407-420. [2] Zaiss M, et al. *NMR Biomed*. **2018**, 31(4):e3879. [3] Feng L, et al. *Magn Reson Med*. **2013**, 70 (1), 64-74. [4] Chen L, et al. *Magn Reson Med*. **2019**; 81(1), 69-78




Fig. 1. images reconstructed using (a) 5 shots, (b) 2 shots and (c) a single shot.




Fig. 2. Z-spectra of an ROI from images reconstructed using 5, 2 and 1 shot(s).




Fig. 3. Top: NOE maps calculated from images reconstructed using 5, 2 and 1 shot(s). Bottom: corresponding Z spectra and the PLOF fittings for a white matter ROI.